

Lipase Inhibitor Investigation

Lipase Inhibitor Investigation

	l able of Content	
1	 Introduction Overview Teaching Plan & Key Features Important Notes 	Page 1-2 1 1-2 2
2	 Instructional Materials Student Worksheet 1 Student Worksheet 2 Student Samples 1 (Worksheet 2) Teacher Notes 1 Laboratory Manual Teacher Notes 	3-19 3-5 6-8 9-11 12-13 14 15-19
3	 Supplementary Resources Possible Modifications Technician Notes References 	20-21 20 21 21



Notes for teachers

• Scan the QR code to get the electronic files.

• Teachers are strongly encouraged to adapt and modify these resources as necessary.





Lipase Inhibitor Investigation

Overview

- The *Lipase Inhibitor Investigation* is about the search for an anti-obesity agent. Students investigate the inhibitory effects of different types of bitter melon seed extracts on lipase activity using the milk–pH indicator system (Royal Society of Biology Nuffield Foundation, 2019).
- Students are given the opportunity to design and carry out experiments in which they set up controls, consider the need for replicates, and identify limitations of using visual inspection to determine the end point of a reaction and an *in vitro* system to study the effects of the seed extracts on enzyme activity *in vivo*.

Teaching Plan & Key Features

Prerequisite knowledge (scientific ideas)

- Food substances and energy requirement in humans
- Digestion and absorption of fats in humans

Prerequisite manipulative skills

• Using an autopipette to transfer a small volume of solution

Lesson	Lesson sequence	Duration (mins)	Resources					
Stage O Preparing for the investigation								
• It is situated in an authentic context related to the search for anti-obesity drugs (Contextualisation).								
Students rea	ad information about the background of the investigation (<i>Reading</i>	g Materials).						
Before	• The teacher distributes <i>Worksheet 1</i> for students to complet	e at home so	Worksheet 1					
Lesson I	that they can be familiar with the background of the investi-	gation.						
1	• The teacher discusses the investigation context with	40	Worksheet 2					
	students.							
	• The teacher provides feedback on students' responses in							
	Worksheet 1.							
	• The teacher distributes <i>Worksheet 2</i> for students to							
Before	• The teacher distributes <i>Worksheet 2</i> for students to complete at home. <i>Worksheet 2</i>							
Lesson 2								
Stage 2 Design	Stage 2 Designing the investigation							
 Students int 	eract with a virtual laboratory to familiarise themselves with the n	materials and a	pparatuses they					
would use i	would use in the investigation (Virtual Laboratory).							
 Students use 	• Students use a template to design their own experimental set-ups (<i>Investigation Planning Template</i>).							
 Students has 	• Students have the chance to evaluate their own and their peers' experimental set-ups (Self & Peer							
Evaluation)								
2	• Teacher provides feedback on students' experimental	40	Student Samples 1					
	designs in <i>Worksheet 2</i> .							
3	• The teacher discusses with the students some questions	40	Teacher Notes 1					
	related to the experimental design.							
	• The teacher provides students with laboratory manual for							
	preparation at home.							

 Stage S Carrying out the investigation Students use microscale instrumentation that reduces the time of the experiments (<i>Microscale Instrumentation</i>). Students collect more complex data sets by setting up replicates (Complex Data Set). Students use camera to collect data (<i>Digital Tool</i>). 								
4	 4 Teacher asks questions to help students connect their lab experience and related ideas/scientific inquiry skills. • Students carry out the investigation. 							
 Stage 9 Explaining and evaluating data Students use <i>Google Sheet</i> for data recording and manipulation (<i>Digital Tool</i>). Students use data to identify seed extracts with the highest inhibitory effect on lipase activity. Students considering the limitations of using an <i>in vitro</i> system to study the effects of the seed extracts on enzyme activity <i>in vivo</i>. 								
Before Lesson 5• Students complete data reporting and analysis at home.Teacher Notes 2• Teacher collects and marks student responses.								
5	• Teacher provides feedback on students' performance related to data reporting and analysis.	40	Teacher Notes 2					

Important Notes

- Students are *not* required to learn the detailed mechanism of enzyme inhibition. Rather, they are expected to use data to support their claims about the inhibitory effect.
- Students should avoid skin contact with the solutions and quickly rinse any splashes of lipase solution or sodium carbonate from their skin.



Instructional Materials

Stage ① Preparing for the investigation

Student Worksheet 1



Notes for teachers

- The teachers can distribute *Worksheet 1* and ask students to read the background information related to the investigation at home.
- Students' responses can be collected using a *Google Form*.

<u>Task 1</u>

- Read the following information and source materials in the *Data File*.
- Answer the questions that follow.

Obesity is a major risk factor of cardiovascular diseases, musculoskeletal disorders, and some cancers. According to the World Health Organization, approximately 650 million adults were obese, and more than 1.9 billion were overweight in 2016.

Orlistat is a drug approved by the United States Food and Drug Administration for the long-term treatment of obesity; it inhibits lipase activity in the alimentary canal. Orlistat reduces the absorption of dietary fat in the human body. However, it may cause side effects such as gastrointestinal discomfort.

Scientists are now searching for natural alternatives as anti-obesity drugs. Read the information in the *Data File* to familiarise yourself with the investigation background.



Scan the QR code to get a copy of the Google Form.



Teachers can diagnose students' difficulties in understanding the relevant content and methods for measuring lipase activity, and then provide feedback before students design the experiments.

<u>Data File</u>

Your biology teacher asks you to read the following source materials to prepare you to design a scientific investigation related to lipase activity:

Source 1:



Source 2:

Student scientists in Hong Kong discover anti-obesity agents present in bitter melon seed extract

Some plant tissues contain lipase inhibitors such as polyphenols and saponins.

A team of secondary school student scientists screened more than 60 plant samples to identify natural inhibitors that effectively reduce pancreatic lipase activity *in vitro*. They performed a simple experiment using whole milk and an alkaline solution containing a pH indicator (which is blue at alkaline pH values and yellow at acidic pH values). They first prepared different seed extracts by grinding the seeds with a pestle in a mortar using water and a

spoonful of sand.

Scan the QR code to access the material used by the student scientists.

The plant extract was then incubated with alkaline pancreatic lipase containing the pH indicator for 5 minutes. After adding whole milk to initiate fat digestion, the time taken for the

reaction mixture to change colour from blue to yellow was recorded. The recorded data were then used to determine the activity of the pancreatic lipase preincubated with different plant extracts.

Scan the QR code to watch their investigation.





The student scientists' findings revealed that bitter melon (*Momordica charantia*) seed extracts contain pancreatic lipase inhibitors. Their findings show potential for addressing the global obesity problem.

Answer the questions below *after* reading the source materials:

- (a) Explain why inhibiting the lipase activity in the alimentary canal can help reduce body weight.
- (b) In which part of the alimentary canal can you find pancreatic lipase? Explain the conditions that favour the pancreatic lipase activity in this part.
- (c) Whole milk contains triglycerides. Write a word equation to show the actions of pancreatic lipase on triglycerides in whole milk.
- (d) Explain why whole milk containing an alkaline solution and the pH indicator described in *Source 2* would turn from blue to yellow after the addition of pancreatic lipase.
- (e) How is the time taken for the alkaline solution to turn from blue to yellow related to the rate of lipase activity?
- (f) After reading the source material, propose *one* investigation question related to the material you have read.

Stage 2 Designing the investigation

Student Worksheet 2

<u>Task 2</u>

Answer the questions that follow. •

Scenario

Bitter melon (Momordica charantia) belongs to the Cucurbitaceae family. Different varieties of bitter melons have different shapes and bitterness. The Cucurbitaceae family is composed of different types of melons.

Your biology teacher asks you to design an investigation to compare the inhibitory effect of the seed extracts of three different types of melons within the Cucurbitaceae family on pancreatic lipase. The goal is to identify the seed extract sample with the highest inhibitory effect on pancreatic lipase activity.

You received the following materials:

Alkaline solution containing a pH	Glass vials	Test sample 1 [Bitter melon 1
indicator (blue under alkaline pHs		(Momordica charantia) seed
and yellow under acidic pHs)		extract]
Pancreatic lipase	Timer	Test sample 2 [Bitter melon 2
		(Momordica charantia) seed
		extract]
Orlistat (a drug that inhibits	Distilled water	Test sample 3 [Angled Luffa
pancreatic lipase)		(Luffa acutangular) seed
		extract]
Whole milk	Tablet (to be used as a camera)	Tablet stand
Autopipette	Autopipette tips	



Scan the QR code to view the materials.



The virtual laboratory provides students with opportunities to get familiar with materials used in the investigation.

(a) Briefly describe how you would use the materials to design an investigation to achieve the aim. You can also draw your experimental design.

(For this purpose, the Investigation Planning Template may be helpful.)



Brief explanation of my design:



Scan the QR code to get a copy of the *Investigation Planning Template*.





| 2.5-8

Student Samples 1 (Worksheet 2)

Notes for teachers

- After collecting students' designs, teachers can select student drawings (anonymised) for discussion.
- The following shows three samples with varying sophistication in responses. Some guiding questions can be included to facilitate students' evaluation of experimental designs.

Examples of students' experimental designs

Possible guiding questions

- Which design(s) accurately represent the independent variable and provide methods for manipulating it? Why does your group think so?
- Which design(s) accurately represent the dependent variable and specify the parameters for measuring it? Why does your group think so?
- Which design(s) demonstrates the correct sequence for adding the chemicals? Why does your group think so?
- Which design(s) incorporate the appropriate control set-up(s), if necessary? Why does your group think so?
- What are the ways to enhance the designs to ensure that the data collected are accurate and reliable?

Design **0**:



Design **Q**:



Design 🛭 :

To investigate the effect of three different helens within the types on the poncreatic lipage Cucurbitalece family 04 activities. Offistet Pistille} panerectic Wetos Mixed Julition Mirel Solution whole will rixed solution whole will , alkeline sulltion. alkeline sulltion, (pencreatic lipere pencientia lipere alkeline Julitine Note: Set a timer for each setup for one minute once the solution in the micropigette tip is added to the mixed solution invide the glass vial. Test Sample 3 Test Sample 2 Sample Test Test sample 3 Test sample 2 (Bitter helon 2) Test Sample 1 (Aude) (Independent Variable) (Bitter nelin 1) (Independent Variable) (Integentate Variable) Mixed Silution niked Julution whole milk, alkeline siletion Mixed sulation while (whole milk, alkoline solution, pencreetic lippoe al feline solution, paherentic lipare parcreatic lipse Note: Set a timer for each setup for one minute the test Sample 044 in the micompipette tip is added to the mixed solution inside the glass vial. Contral variables: Thitled poncreative lipse and concentration of test sample. This is to easure the accuracy of the experiment. enter the solution through air during experiment. Barteria may

Notes for teachers

• Teachers can capture and represent student thinking using public displays (e.g. whiteboards) and then work with students to explore their divergent thinking.

Group	Laura	Zee	Nicol	e Kell	y Cand	Cindy	Iris	Roiko	(Lily	Sonia	
RI	1.3	3	2,3	1,3	1.3	3	1,3			1.3	
27	×	- 1	1.2	2	2	2	1	2	1,2	1,2	
0 1.											
Q 5.											
Q. 4.											

Teacher Notes 1



Notes for teachers

- After receiving feedback on their experimental designs, the following shows questions that teachers may use to guide students in thinking about and assessing the scientific inquiry skills related to their experimental designs.
- Student work samples are shown below to illustrate possible student thinking to some questions.

<u>Task 3</u>

Possible questions

1. The following shows two methods to measure the dependent variable:

Method A:	Record the colour of the alkaline solution containing the pH indicator in glass vials, both <i>with</i> and <i>without</i> the test samples, after 10 min.
Method B:	Measure the time it takes for the colour of the alkaline solution containing the pH indicator to change (i.e. reach the end point) in the glass vials, both <i>with</i> and <i>without</i> the test samples.

Your teacher suggests that you should use *Method B*.

- (a) What is the limitation of using *Method A* to compare the inhibitory effects of different types of seed extract on pancreatic lipase activity?
- (b) Explain how the inhibitory effects of different types of seed extract on pancreatic lipase can be compared using *Method B*.
- 2. Jeffrey proposes two methods for setting up a control to compare the degree of inhibition on alkaline lipase activity:
 - *Set-up A*: Replacing the test samples with orlistat.
 - *Set-up B*: Replacing the alkaline lipase with boiled alkaline lipase in the glass vials containing the test samples.

Which set-up, A or B, enables a more accurate determination of the degree of inhibition? Explain your answer.

(Put a ' \checkmark ' in the appropriate box.)

I will choose \Box Set-up A \Box Set-up B

The reasons:



Notes for teachers

- Q.1 assesses students' ability to connect the methods of measurement to the dependent variable and the limitations related to the measurement method.
- Q.2 assesses students' ability to set up the control and explain its function.

The following are some examples of students' responses to Q.1(a):

<u>Sample 1</u>

UZBDGDED (1) Your teacher suggests that you should use *Method B*. What is the limitation of using Method A to compare the inhibitory effect of different types of seed extract on pancreatic lipase activity? The colour present is based on personal tory effect of different to compare the on pancreatic seed extract lipase activity Sample 2 UDBDGDED Your teacher suggests that you should use Method B. (1) What is the limitation of using Method A to compare the inhibitory effect of different types of seed extract on pancreatic lipase activity? We cannot compore the inhibitory effect of different types of reed extract on pancreatic lipase activity if there is no colour change in pH indicator within to minutes by using Method A.

<u>Sample 3</u>

(1) Your teacher suggests that you should use Method B. What is the limitation of using Method A to compare the inhibitory effect of different types of seed extract on pancreatic lipase activity 26 Investigate lipase activity for method A1 we cannot compare the colour Intensity of set-ups IF LO minutes, enough for lipase to digest all lipids, which will give same or similar colour threesity of all set-ups. We also cannot compare the colour Intensity of all sample change from blue to yellow if no set-ups have colour change after to minutes. Thus, we cannot compare the deaves of unitorities of lipase autivity.



About the samples

- Sample 1 incorrectly cited the limitation arising from the subjectivity of colour judgment, a limitation inherent in both methods.
- Sample 2 correctly identified the limitation but lacked details whereas Sample 3 provided a detailed explanation, such as the absence of colour difference because all lipids were digested within the specified time frame.

Laboratory Manual



<u>Task 4</u>

Read the following procedures to carry out the investigation.

Procedure

Determining the inhibitory effects of seed extracts

- 1. Use your mobile phone/tablet to start recording a video.
- 2. Label 15 vials (A1–3 to E1–3).
- 3. Add the seed extracts/orlistat/distilled water, alkaline solution containing the pH indicator, and pancreatic lipase into the vials, according to the following table:

Vial	Sample (mL)		Alkaline solution containing the pH indicator (mL)	Pancreatic lipase (mL)
Α	Seed extract 1	1	2	2
В	Seed extract 2	1	2	2
С	Seed extract 3	1	2	2
D	Orlistat	1	2	2
Е	Water	1	2	2

- 4. Incubate the vials at room temperature for 5 minutes.
- 5. Add 3 mL of whole milk to the vials, and shake the vials well.
- 6. Start the timer.
- 7. Shake the vials occasionally.
- 8. Repeat *Steps 4–8* two more times.
- 9. Stop the video recording when the colour of the solution in all of the vials turns yellow from blue.
- 10. Fill in the data in the *Google Sheet*.



Scan the QR code to get a copy of the Google Sheet.



The *Google Sheet* helps students process and visualise the data they collected.

Teacher Notes 2



Notes for teachers

- The following are possible questions that teachers can use to guide students in thinking about or assessing their scientific inquiry skills related to data analysis and interpretation.
- Student work samples are shown below to illustrate possible student thinking to some questions.

<u>Task 5</u>

Possible questions

1. Plot a graph to show the effect of the seed extracts on pancreatic lipase activity.

- Consider the following questions when plotting the graph:
 - Which type of graph (bar graph, line graph, pie chart, etc.) would you choose? Why?
 - Which axis (x-axis/y-axis) should contain the independent variable?
 - Which axis (x-axis/y-axis) should contain the dependent variable?
 - What would be a suitable title for your graph?

Some reminders are added to guide students in constructing graphical representations appropriately.

- 2. Boris found that an outlier was present in a replicate of one of the test samples. Suggest *one* possible reason for why this occurred.
- 3. Your classmate claims that the seed extract that showed the highest inhibitory effect on pancreatic lipase in this investigation should be used as an anti-obesity drug.

Discuss whether you agree with this claim. (Put a ' \checkmark ' in the appropriate box.)

□ Agree

- Disagree
- 4. Suggest one *new* investigation that needs to be conducted before the seed extract(s) can be used as an anti-obesity drug.



Notes for teachers

- Q.1 assesses students' ability to construct appropriate graphic representations.
- Q.2 assesses students' ability to propose reasons to explain the occurrence of an outliner.
- Q.3 assesses students' ability to identify the limitations of the generalisability of the results from an *in vitro* to an *in vivo* system.
- Q.4 assesses students' ability to generate a new investigation question that extends the present investigation.

The following are some examples of students' responses to Q.1:

<u>Sample 1</u>





About the samples:

- Sample 1 mistakenly used a line graph for data presentation while Sample 2 correctly used a bar chart to present the data. A bar chart is appropriate as the independent variable is the type of seed extract, which is a categorical variable.
- Sample 1 also has a proper title and labelling of the x and y axes.

The following are some examples of students' responses to Q.2:

<u>Sample 1</u>

(2) Boris found that an outlier was present in one of the replicates of one of E2: U B B G B E B the test samples. Suggest one possible reason for why this occurred.
Mult versult in a great diffume in Mult versult in a great diffume in Mult versult in a great diffume in the colour change is difficult to observe. We can't a great the data? exactly tell when the colour chang. The time
(2) Boris found that an outlier was present in one of the replicates of one of E2: U B B G B E B B G B E

The sample used does not have equal amount of inhibitary effect within the same sample used. The part of sample used the sit up due to have less inhibitar than other taken has much more or much less inhibitar than other part of the same sample.

<u>Sample 3</u>

(2) Boris found that an outlier was present in one of the replicates of one of $E2: U \square B \square G \square E \square$ the test samples. Suggest *one* possible reason for why this occurred.

There might be human error occurs. The person might have made mistake in adding a different amount of milk into jneans a different our out of Imidis added, the the sample, which requires a longer or shorter period of time for the breakdown of lipid and the colour change.



About the samples

- Sample 1 proposed a reason that was not sufficiently convincing in explaining the occurrence of the outlier.
- Sample 2 proposed a plausible reason for the occurrence of the outliers but lacked a detailed explanation.
- Sample 3 not only provided a plausible reason but also offered a more thorough explanation.

The following are some examples of students' responses to Q.3:

<u>Sample 1</u>

Discuss whether you agree with this claim. (Put a "1" into the appropriate box.) EII: UDBDGDED □ Agree M Disagree that Invertigation Carrier out man · Not IVO

<u>Sample 2</u>

Discuss whether you agree with this claim. (Put a " \checkmark " into the appropriate box.) E11: U \Box B \Box G \Box E \Box

□ Agree Disagree a experiment untside we donot donot chongh 15 as e Ffeet know its will donot) this bud y the substance in body each as (Stu) temperature body , different situation in body, the effect cannot reflect in our body. ouid totally different STUMACh extract May Do



About the samples

- Both samples correctly disagreed with the claim. Both samples identified the limitations of generalising results produced in an *in vitro* system to the *in vivo* conditions.
- Sample 2 further provided an explanation for the difference between the *in vivo* and *in vitro* conditions.

The following are some examples of students' responses to Q.4:

<u>Sample 1</u>

(i) Suggest *one* new investigation that needs to be conducted before the seed G12: U □ B □ G □ extract(s) in (g) can be used as an anti-obesity drug.

Investigate whether have side effect on human's U J TENT

<u>Sample 2</u>

(i) Suggest *one* new investigation that needs to be conducted before the seed extract(s) in (g) can be used as an anti-obesity drug.

Carry our the whole investigation in 1200 (e.g. white rat) or at a set temperature of 37°C

<u>Sample 3</u>

(i) Suggest *one* new investigation that needs to be conducted before the seed $G12: U \square B \square G \square$ extract(s) in (g) can be used as an anti-obesity drug.

To investigate whether the extracts will inhibit other like anylase like anylase enzymes, present in pancreatic juice.



About the samples

• All the samples were able to generate a new investigation question that expanded upon the current investigation. However, the question posed by Sample 1 was somewhat vague.

Supplementary Resources

Possible Modifications

1. Using immobilised lipase beads to investigate lipase activity

- Lipase can be immobilised using 3% sodium alginate solution. Immobilised lipase beads can be used to study lipase activity using a milk-pH indicator system.
- The following shows the sample results of the investigation that examined the effect of increasing the number of lipase beads on the digestion of milk.
- See Chan et al. (2024) for procedures of how to make immobilised enzyme beads.



Technician Notes

1. Materials for Task 4

Preparation of seed extract

- 1. Weigh 2 g of seeds using an electronic balance.
- 2. Place the seeds in a mortar and pestle.
- 3. Add a spoonful of sand.
- 4. Add 10 mL of distilled water.
- 5. Grind the seeds into powder.
- 6. Filter the extract/Centrifuge the extract at top speed (13, 500 rpm) to obtain the supernatant.

(A grinder can be used to grind the seeds.)

Chemicals to be prepared

- Alkaline solution with a pH indicator (a master mix comprising 100 mL of 2% sodium carbonate [2 g of sodium carbonate in 100 mL distilled water] and 200 mL of 0.04% bromothymol blue [0.1 g of bromothymol blue in 16 mL of 0.01 M sodium hydroxide, with the volume made up to 250 mL with distilled water])
- 5% porcine pancreatic lipase (0.5 g in 10 mL distilled water)
- 10 mg/mL orlistat (120 mg tablet dissolved in 1 mL of absolute ethanol, with the volume made up to 12 mL with distilled water)

Mat	terials for each group		
•	Whole milk (>35 mL)	• Vials X 15	• Tablet stand
•	Distilled water (>3 mL)	• Timer	Tablet/mobile phone
•	Seed extract 1 to 3 (>3 mL)	• Orlistat (>3 mL)	Rubbish bin
•	Autopipette (P-1000)	• Autopipette tips (P-1000)	• Labels
•	Pen	• Alkaline solution with a pH	Pancreatic lipase
		indicator (>30 mL)	solution (>30 mL)



References

- Chan, K. K. H., Ho, D. T. S., & Lau, D. S. P. (2024). Using amylase beads to investigate factors affecting enzyme activity. *The American Biology Teacher*, *86*(3), 153–160.
- Chan, P. C., & Chan, K. K. H. (2023). Inquiry on a potential anti-obesity agent: Investigating pancreatic lipase inhibitors in seed extracts *The American Biology Teacher*, *85*(5):265–269
- Royal Society of Biology Nuffield Foundation. (2019). Investigating effect of temperature on the activity of lipase. https://practicalbiology.org/bio-molecules/factors-affecting-enzyme-activity/investigating-effect-of-temperature-on-the-activity-of-lipase

This project is supported by Quality Education Fund (Project No. 2019/0283)

The Trustee of the Quality Education Fund is the owner of the copyright of this product. Any reproduction of this product for commercial purposes is strictly prohibited unless prior written consent has been obtained from the Trustee of the Quality Education Fund. 此計劃由優質教育基金贊助(計劃編號 2019/0283) 此計劃產品版權屬優質教育基金擁有,未經許可, 不得翻印以作商業用途。

侵質教育基金





